



150250

Tork Couch Roll Advanced

Protect your patient couch with the disposable Tork Couch Roll Advanced. The perforated couch roll is ideal to maintain a high level of hygiene in a healthcare environment. Rolls can be used in Tork Couch Roll Dispenser.

- Perforation for easy tear off and consumption control
- Strong and easy to roll, for optimal protection
- Better resistance and higher comfort due to 2-plies



Tork Couch Roll Advanced

Product data

System C1

Roll width 50 cm

Sheet legth 33 cm

Number of Sheets

165

Quality Advanced

Roll Diameter 12 cm

Core Inner Diameter

5.1 cm

Embossing Yes Roll length 54.5 m

Number of Sheets

Ply 2

Colour White

Product certificates





The mark of responsible forestry

Compatible products





129184 Tork Couch Roll Dispenser

129183 Tork Couch Roll Dispenser

Shipping data

| | Consumer units (CON) | Transport unit (TRP) | Pallet (PAL) |
|------------------|----------------------|-------------------------|-----------------|
| EAN | 7310791230351 | 7310791230344 | 7322540904192 |
| Packing material | Plastic | Plastic | |
| Pieces | 1 | 9 (9 CON) | 162 (18 TRP) |
| Gross weight | 1025.3 g | 9.3 kg | 167.54 kg |
| Height | 500 mm | 500 mm | 1161 mm |
| Length | 120 mm | 360 mm | 1200 mm |
| Width | 120 mm | 360 mm | 1000 mm |
| Net weight | 980.1 g | 8.8 kg | 158.78 kg |
| Volume | 7.2 dm3 | 64.8 dm3 | 1.17 m3 |